

# QYNAPSE


**QYPREDICT**<sup>®</sup> PROGNOSTIC MODEL ENRICHES FOR FASTER DECLINERS IN MILD COGNITIVE IMPAIRMENT

**Jorge Samper-Gonzalez**, Elizabeth Gordon, Enrica Cavedo, Clarisse Longo dos Santos, Adam J. Schwarz

Qynapse

**PARIS, FRANCE & BOSTON, USA**

AD/PD 2022 Presentation – March 18<sup>th</sup> 2022







# AD/PD™ 2022

ADVANCES IN SCIENCE & THERAPY

International Conference on  
Alzheimer's and Parkinson's Diseases  
and related neurological disorders

March 15 - 20, 2022 | Barcelona, Spain **+ HYBRID**



## DISCLOSURES

	Nothing to disclose
<b>X</b>	Yes, please specify

Company Name	Honoraria / Expense	Consulting / Advisory Board	Funded Research	Royalties / Patent	Stock Options	Ownership / Equity Position	Employee	Other (Please specify)
Qynapse							<b>X</b>	

# Background

## Key Challenges:

- Many patients enrolled in clinical trials for Alzheimer's Disease (AD) interventions do not progress clinically over the study period, reducing the power to detect positive treatment effects.
- The increased heterogeneity and slower decline in mild cognitive impairment (MCI) pose even greater challenges but this stage is a better window to target intervention.
- The suboptimal selection of patients has been a key contributor to the reduced success rate of disease-modifying trials and improvements in selection strategy are urgently needed to better power clinical trials.

## New Approaches:

- Recent advances in AI predictive modeling, such as the QyPredict® algorithm, are promising tools to improve the selection of patient populations more likely to clinically progress during the timeframe of an AD clinical trial.

## Study objectives

1. To evaluate the prognostic value of QyPredict® in all-comer mild cognitive impairment (MCI) populations as well as amyloid-positive ( $A\beta+$ ) and APOE- $\epsilon$ 4 + sub-populations.
2. To evaluate the benefit of using QyPredict® to refine patient selection in terms of clinical trial success probability in a simulation study.




# INTRODUCTION

## QYPREDICT®

### What is QyPredict®?


- QyPredict® is a prognostic model currently developed for use in mild cognitive impairment (MCI) populations
- QyPredict® takes different **baseline** inputs : structural MRI outputs from QyScore®, demographic and clinical data, genetic and biological disease markers





### What information does QyPredict<sup>®</sup> provide?

- For each individual, QyPredict<sup>®</sup> generates a score (between 0-1) representing a probability of having a specific outcome (e.g. an increase of at least 0.5 points of CDR-SOB score over 24 months)
- Allows for a personalized medicine approach by providing prediction on an individual patient level





## Participants

Study participants were from **ADNI** with the following inclusion criteria:

- Age (55 – 85 years of age)
- MCI diagnosis
- MMSE between 24 and 30
- CDR = 0.5 at baseline
- Available amyloid status



## Analyses

- A QyPredict® probability of decline was calculated for each of the 519 individuals modeled **over 24 months.**
- The performance of QyPredict® to accurately model real decline in CDR-SOB was evaluated at several **increasing probability of decline cut-offs of 0.1, 0.2, 0.3, 0.4 and 0.5.**
- For the total cohort, Amyloid+ and APOE-ε4+ populations, the mean, standard deviation and Cohen's d for CDR-SOB change was calculated.
- Sensitivity, specificity, positive predictive value (PPV) and power calculations were computed to assess predictive performance.
- Finally, a clinical trial simulation was run to further investigate the utility of QyPredict® to improve the probability of clinical trial success.



# Results

Demographics and CDR-SOB for the Full and filtered QyPredict® Cohorts


- Age increased with increasing QyPredict® threshold ( $p < 0.001$  only for QyPredict® > 0.4)
- Baseline CDR-SOB increases but it is only significant when comparing the full and QyPredict® > 0.5 groups
- Change in CDR-SOB significantly increased based on the baseline QyPredict® values
- Proportion of Amyloid positive individuals increases above QyPredict® of 0.2

Cohort	N	Age	Sex (M/F)	CDR-SOB at baseline	CDR-SOB at 12 months	CDR-SOB at 24 months	Amyloid positive N (% total)	APOE-ε4 status (+ / -) (% +)
Full cohort	519	71.8 (±7.1)	303 / 216	1.5 (±0.90)	1.8 (±1.38)	2.3 (±2.11)	319 (61%)	263 / 256 (51%)
QyPredict® > 0.1	467	72.2 (±7.0)	287 / 180	1.6 (±0.91)	1.9 (±1.39)	2.5 (±2.14)	300 (64%)	243 / 224 (52%)
QyPredict® > 0.2	418	72.6 (±7.0)	266 / 152	1.6 (±0.90)	2.0 (±1.4) *	2.7 (±2.17) *	283 (67%) *	223 / 195 (53%)
QyPredict® > 0.3	353	73.2 (±6.8)	224 / 129	1.6 (±0.93)	2.2 (±1.41) *	3.0 (±2.21) *	257 (72%) *	195 / 158 (55%)
QyPredict® > 0.4	313	73.4 (±6.6) *	199 / 114	1.7 (±0.95)	2.3 (±1.42) *	3.1 (±2.25) *	234 (74%) *	179 / 134 (57%)
QyPredict® > 0.5	251	73.9 (±6.6) *	154 / 97	1.7 (±0.97) *	2.4 (±1.44) *	3.4 (±2.30) *	197 (78%) *	148 / 103 (59%)




# Results: Stable versus Decliners

- QyPredict® performed well at predicting those that would decline and those that would remain stable over 24 months in all Cohorts.
- Stable defined as a QyPredict® score < 0.5 and Decliner as QyPredict® > 0.5




**251** subjects were predicted to decline.  
**190** actually declined, for a **PPV\*** of **0.76**,  
**sensitivity** of **0.7** and **specificity** of **0.75**

Sample size	Un-Enriched Cohort	Enriched Cohort
30% treatment effect	304	102
50% treatment effect	111	38



**197** subjects were predicted to decline.  
**164** actually declined for a **PPV** of **0.83**, **sensitivity**  
of **0.77** and **specificity** of **0.70**

Sample size	Un-Enriched Cohort	Enriched Cohort
30% treatment effect	144	77
50% treatment effect	53	29



**148** subjects were predicted to decline.  
**122** actually declined for a **PPV** of **0.82**,  
**sensitivity** of **0.75** and **specificity** of **0.74**

Sample size	Un-Enriched Cohort	Enriched Cohort
30% treatment effect	184	79
50% treatment effect	68	30

\* PPV = positive predictive value



# QYPREDICT<sup>®</sup> and Clinical Trials Simulation

QYNAPSE

QYPREDICT<sup>®</sup>  
incorporated into  
the inclusion criteria  
at screening

Original database ADNI  
1987 patients

Inclusion  
criteria

- Age: [55-85]
- MMSE: [24-30]
- CDR: 0.5
- QYPREDICT<sup>®</sup>
- Not used for model training

Eligible patients  
143 patients

Bootstrap with replacement

Selected patients  
1000 patients

Control arm  
500 patients

Treatment arm  
500 patients

Statistical analysis


Success /  
Failure

Repeated for  
1000 simulations

Treatment effect:  
30% reduction on  
CDR-SOB decline

# Clinical Trial Simulation Results

QYPREDICT® incorporated into the inclusion criteria at screening



n = total patients screened to reach 1000 patients enrolment



## SUMMARY AND CONCLUSIONS


- Using baseline neuroimaging and demographic information, QYPREDICT® was able to accurately model the likelihood an individual patient would decline over 24 months in all-comers, Amyloid positive and APOE positive populations, based on change in CDR-SOB.
- Sensitivity, specificity and positive predictive value were high across the different populations (all > 0.70).
- Enriching using QYPREDICT® substantially reduced sample sizes required to detect a treatment effect.

## SUMMARY AND CONCLUSIONS

- QYPREDICT® shows promise in improving trial selection towards decliners for increased trial success probability with a single upfront screening cost.
- The use of QYPREDICT® score as part of the inclusion criteria in our clinical trials simulation **significantly** improved the probability of trial success, while increasing screening failure rates due to excluding those who would be less likely to clinically progress.
- These results support the promising potential to improve design and power of AD clinical trials, and the likelihood of **detecting positive treatment effects** and **achieving trial success**.



THANK YOU FOR YOUR ATTENTION AND  
I WELCOME ANY QUESTIONS



FOR MORE INFORMATION:

Jorge Samper-Gonzalez  
Data Scientist  
Qynapse  
Email: [jsamper@qynapse.com](mailto:jsamper@qynapse.com)





PEACE OF MIND