

# Use of predictive algorithms for the selection of patients in clinical trials: an enrichment strategies comparison



Antoine Movschin<sup>1\*</sup>, Clarisse Longo dos Santos<sup>1</sup>, Achille Mascia<sup>1</sup>, Jorge Samper-González<sup>1</sup>, Urielle Thoprakarn<sup>1</sup>, Philippe Tran<sup>1,2</sup>, Jean-Baptiste Martini<sup>1</sup>, Enrica Cavedo<sup>1</sup>, Alzheimer's Disease Neuroimaging Initiative

<sup>1</sup> Qynapse SAS, Paris France (<https://www.qynapse.com/>) <sup>2</sup> Equipe-projet ARAMIS, ICM, CNRS UMR 7225, Inserm U1117, Sorbonne Université UMR\_S 1127, Centre Inria de Paris, Groupe Hospitalier Pitié-Salpêtrière Charles Foix, Faculté de Médecine Sorbonne Université, Paris, France

\* amovschin@qynapse.com

## Background

---

The majority of disease-modifying trials (DMT) for Alzheimer's disease (AD) drugs have failed since 2003 (Yiannopoulou KG et al., Biomedicines. 2019)

---

Suboptimal patient cohort selection is one of the main causes for high trial failure rates in AD

---

Recent advances in artificial intelligence (AI) could help select Mild Cognitive Impairment (MCI) patients who are more likely to decline cognitively

## Objectives


---

To compare various recruitment strategies for patient enrichment in terms of trial success probability and screening failure rate

---

To evaluate the benefit of using a predictive tool such as QyPredict<sup>®</sup> to refine patient selection


# Methods



**decliner:** MCI patient with an increase > 0.5 in CDRSOB at 24 months  
**non-decliner:** MCI patient with no increase > 0.5 in CDRSOB at 24 months  
**CDRSOB:** Clinical Dementia Rating – Sum Of Boxes  
**CDRSOB APC:** CDRSOB Annual Percent Change  
**MCI:** Mild Cognitive Impairment

**MMSE:** Mini-Mental State Exam  
**CDR:** Clinical Dementia Rating  
**QyPredict®:** software that predicts whether an MCI subject is going to have a strictly positive increase in CDRSOB at 24 months

# Results



## Conclusion

We presented a framework to evaluate the success probability and screening failure rate of various patient selection strategies for a clinical trial.

Such a framework could also be used to derive the cost and the duration of a trial, depending on the patient selection strategy.

The use of a predictive solution such as QyPredict® to enrich the patient selection strategy increases the probability of success of the trial.

A trade-off is to be made, taking into account the trial's objectives in terms of success likelihood, cost and duration.