

# FOUR-YEAR COGNITIVE DECLINE IN MILD COGNITIVE IMPAIRMENT PATIENTS WITH VARYING BASELINE NEUROANATOMICAL ATROPHY PROFILES: A MEMENTO AND QYSCORE® STUDY.

E. Gordon<sup>1</sup>, A.J. Schwarz<sup>2</sup>, E. Cavedo<sup>1</sup>, C. Longo dos Santos<sup>1</sup>, J. Samper-Gonzales<sup>1</sup>, and the MEMENTO Study Group


QYNAPSE

<sup>1</sup>Qynapse, Paris, France <sup>2</sup>Qynapse, Boston, MA. - email: egordon@qynapse.com

## BACKGROUND

Three patterns of brain atrophy have been consistently identified in cases of Alzheimer's dementia<sup>1</sup> and may provide important prognostic information about future cognitive decline in mild cognitive impairment (MCI)<sup>2,3</sup>

- **limbic-predominant (LP)**: atrophy of the hippocampus but relatively spared neocortex
- **hippocampal sparing (HS)**: atrophy in the cortex (especially lateral temporal and parietal lobes) but sparing of the hippocampus
- **typical/diffuse (TD)**: atrophy in both hippocampus and cortex



## OBJECTIVES

To assess differential patterns of cognitive decline over 48 months in the MEMENTO population using a volumetric quantitative approach to define atrophy subtypes using the automated output of QyScore® an FDA-approved and CE-marked image analysis software<sup>4</sup>.

## METHODS

### IMAGING DATA

- **N=1153 3DT1 images** from mild cognitive impairment (MCI) participants from the MEMENTO cohort were included (Age=71±8, 688 (59%) female; baseline MMSE = 27.8±2.0, CDR-SB = 0.6±0.7)
- **N=461 (40%)** had amyloid status from CSF or PET, with **145 amyloid-positive (A+)**.
- Neurodegeneration-positive (**N+**) atrophy subtypes (**LP, HS or TD**) were defined on individual scans using z-scores from parietal cortex and hippocampus volumes referenced to 1290 cognitively normal individuals using **QyScore® v1.10**
- Volumes were corrected by **intracranial volume (ICV)**.

Atrophy subtype	N- (no or minimal atrophy)	LP	HS	TD	TOTAL
Total Cohort	1012	37	84	20	1153
Amyloid positive	114	8	21	2	145

### NEUROPSYCHOLOGY DATA AND ANALYSIS

- All **1153** participants underwent cognitive and functional assessments **at baseline** and **48 months**, including mini-mental state examination (**MMSE**), clinical dementia rating scale (**CDR**), Free and Cued Selective Reminding Test (**FCSRT**) and Instrumental Activities of Daily Living (**IADL**) scale.
- Linear regression analysis (correcting for age, sex and baseline performance) investigated the relationship between the baseline neuroanatomical subtype and subsequent 4-year cognitive decline. Bonferroni correction was applied.


## CONCLUSIONS

- This work highlights the utility of quantifying early atrophy profiles in MCI to provide insights into the domains and severity of patient progression over 4 years.


## RESULTS

### NEUROANATOMICAL PROFILE RESULTS

- At **Month 48**, all subtypes demonstrated significant ( $p<0.001$ ) decline across **memory, orientation, and executive function**.
- Severity of decline varied across subscales and subtype and was particularly **marked in A+ LP individuals for memory and orientation**.
- IADL remained stable for the N- and TD subgroups.
- 143 participants converted to **AD [N- = 102 (10%), LP = 18 (49%), HS = 15 (18%), and TD = 8 (40%)** subtypes
- Reflecting these fewer relative conversions, **HS** and **N-** individuals demonstrated the **least decline across all domains** and did not significantly differ on any tests ( $p>0.003$  Bonferroni-corrected).
- **LP atrophy subtype** demonstrated the **greatest decline** across all **individual** cognitive domains as well as global CDR and **IADL** ( $p<0.001$  compared with N- and HS subtypes).
- Similarly, **TD phenotypes** showed significantly decline across multiple cognitive subscales compared with N- and HS subtypes ( $p<0.001$ ) particularly for orientation and global CDR rating.



FCSRT_FREE	Free recall trial of the FCSRT (Memory)
FCSRT_CUED	Cued recall trial of the FCSRT (Memory)
FCSRT_TOTAL	Total recall trial of the FCSRT (Memory)
FCSRT_DELAYED	Delayed recall trial of the FCSRT (Memory)
MMSE_RECALL	MMSE recall section (Memory)
MMSE-OT	MMSE Orientation in time (Orientation)
MMSE-OS	MMSE Orientation in space (Orientation)
MMSE_AC	MMSE Attention and Calculation (Executive function)
MMSE_TOTAL	TOTAL MMSE score (/30) (Global cognition)
CDR-SB	TOTAL CDR-sum of boxes (Global cognition) *
IADL	TOTAL IADL (functional measure of daily living) *



CDR 0 = no dementia,  
CDR 0.5 = questionable dementia,  
CDR 1 = MCI,  
CDR 2 = moderate cognitive impairment and  
CDR 3 = severe cognitive impairment

\* Increasing values for all CDR scores and the IADL indicate greater deficit over time in contrast to declining score on the other neuropsychological assessments

References:  
[1] Ferreira et al. (2020) Neurology 94:436;  
[2] Dufouil et al. (2017) Alz Res Ther 9:67;  
[3] Planche et al. (2021) Alz & Dem 17(4):65;  
[4] Cavedo et al. (2022) Eur Radiol 32(5):2949

**Figure 1.** A) Mean baseline and Month 48 scores on the cognitive and functional assessments, expressed as % of total possible score for each subscale. B) Absolute mean CDR rating score at baseline and Month 48, all presented by baseline atrophy subtype: N- (no/minimal atrophy), LP (limbic predominant), HS (hippocampal sparing), TD (typical/diffuse) A+ (amyloid positive).