

AN AUTOMATED PIPELINE FOR CENTILOID QUANTIFICATION OF AMYLOID- β USING MULTIPLE ^{11}C -PIB-PET AND ^{18}F -PET TRACERS

Elizabeth Gordon¹, Mathilde Borrot¹, Ayoub Gueddou¹, Luca M. Villa¹, Thomas Jubault¹, Nicolas Guizard¹

¹Qynapse, 2-10 rue d'Oradour-sur-Glane, 75015 Paris, France

Email egordon@qynapse.com

BACKGROUND

- Quantitative measures of amyloid- β (A β) pathology using positron emission tomography (PET) imaging are sensitive to identify pathological changes, even at the earliest stages of Alzheimer's disease (AD)
- However, the quantification values vary considerably between tracers and acquisitions, making comparisons across studies and clinical trials findings problematic
- The Centiloid scale aims to standardize these in vivo amyloid quantifications to a 100-point scale, where an average value of zero signifies high certainty of amyloid negativity and 100 identifies average typical AD A β -pathology load¹
- A fully automated Centiloid quantification pipeline supporting multiple available amyloid-PET tracers would be valuable for improving the efficacy and comparability of PET-based analyses across study site

OBJECTIVES

- To develop and validate Qyscore®'s single fully automated Centiloid quantification pipeline for multiple amyloid PET tracers.

METHODS

- Qyscore®'s fully automated pipeline was validated on ^{11}C -PiB-PET and ^{18}F -PET images from the Centiloid project (<https://www.gaain.org/centiloid-project>): 34 young controls [age=31.5 \pm 6.3 years] and 45 AD patients (age=67.5 \pm 10.5 years; CDR= 0.5–1)
- ^{18}F tracers included Flortetapir² (FBP, $n = 46$), Forbetaben³ (FBB, $n = 35$), Flutemetamol⁴ (FTM, $n = 74$) and NAV4694⁵ (NAV, $n = 55$).

- PET/MR image pairs were both co-registered and normalized in the MNI template space (Figure 1).

- The fully automated segmentation from Qyscore®, a CE-marked and FDA-cleared neuroimaging medical device, parcellated the regional masks of the grey matter tissue (target) and of the cerebellum (reference) region (Figure 1)¹

- The standardized uptake value ratio (SUVR) was computed as the ratio of the mean signal in both regions. Correlations of (^{11}C -PiB and ^{18}F) SUVR values with published SUVR data were computed^[2-5].

- Further, correlations between ^{18}F SUVR and paired ^{11}C -PiB SUVR were computed. Correlation coefficients (R^2) > 0.7 were required to consider the Centiloid calibration valid.
- Equations for converting F18-SUVR values to CL were then derived.




Figure 1. Example of Qyscore®'s grey matter composite (target) and cerebellum (reference) masks overlaid onto the MNI transformed PET imaging

RESULTS

Qyscore®'s fully automated quantitative pipeline produced SUVR values well within the bounds defined by the Centiloid method

- SUVR_AD-100 = 2.08 \pm 0.2 and
- SUVR_YC-0 = 1.01 \pm 0.05,
- $R^2 = 0.99$; slope = 1.00; intercept = -0.44).

Qyscore®'s ^{11}C -PiB SUVR correlation coefficients with published values were above 0.99.

Correlation coefficients of Qyscore®'s ^{11}C -PiB SUVR and ^{18}F tracer SUVR's were :

- 0.91 for Flortetapir,
- 0.95 for Forbetaben,
- 0.96 for Flutemetamol,
- 0.99 for NAV4694 (**Figure 2.**)

Equations for converting Qyscore®'s automated SUVR to Centiloid were (Figure 3)

Flortetapir:

$$\text{CL} = 177.79 * \text{FBP_SUVR} - 183.56$$

Forbetaben:

$$\text{CL} = 153.08 * \text{FBB_SUVR} - 152.93$$

Flutemetamol:

$$\text{CL} = 122.39 * \text{FTM_SUVR} - 120.97$$

NAV4694:

$$\text{CL} = 90.20 * \text{NAV_SUVR} - 91.61$$

References: ¹Klunk WE et al. 2015; ²Navitsky M et al. 2018; ³Rowe CC et al. 2017; ⁴Battle MR et al. 2018; ⁵Rowe CC et al. 2016




Figure 2. Linear regression comparing Qyscore®'s automated SUVR for ^{11}C -PiB and the ^{18}F tracers




Figure 3. Qyscore®'s automated SUVR and corresponding Centiloid values for each of the four ^{18}F tracers

CONCLUSIONS

We demonstrate the feasibility and reliability of Qyscore®'s fully automated amyloid PET pipeline for multiple amyloid-PET compounds (PiB and ^{18}F) and transformation to standardized Centiloid quantifications, suitable for implementation in clinical trials.