



**ESR** EUROPEAN SOCIETY  
OF RADIOLOGY

EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*

JULY 13-17, VIENNA & ONLINE

VISIT US AT



MYESR.ORG

# Comparison of the accuracy of quantitative automatic and radiological imaging markers in distinguishing Parkinson Disease and Progressive Supranuclear Palsy

Mambrin F., Cavedo E., Gordon E., Guizard N., Mansueto G., Tinazzi M., Pizzini F.B.

University of Verona – Italy  
Qynapse SAS – Paris (France)



EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*  
JULY 13-17, VIENNA & ONLINE

MYESR.ORG

I do not have any conflict of interest to disclose



EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022


*BUILDING BRIDGES*  
JULY 13-17, VIENNA & ONLINE

MYESR.ORG

# Introduction

## Progressive Supranuclear Palsy (PSP):

- Vertical gaze palsy
- Unsteady gait
- Frequent falls
- Speech and cognitive impairment



# Introduction

Most reliable indexes: M/P ratio, MRPI/MRPI 2.0

MRPI	MRPI 2.0
$MRPI = \frac{P \times MCP}{M \times SCP}$	$MRPI\ 2.0 = MRPI \times \frac{3rdV}{FH}$
P: pons area; MCP: medium cerebellar peduncle width; M: midbrain area; SCP: superior cerebellar peduncle width; 3rdV: 3 <sup>o</sup> ventricle width; FH: frontal horns of lateral ventricles	

Quattrone A, Morelli M, Nigro S, Quattrone A, Vescio B, Arabia G, et al. A new MR imaging index for differentiation of progressive supranuclear palsy-parkinsonism from Parkinson's disease, Parkinsonism and Related Disorders (2018)



EUROPEAN CONGRESS OF RADIOLOGY  
**ECR 2022**

**BUILDING BRIDGES**  
JULY 13-17, VIENNA & ONLINE

MYESR.ORG

# Purpose

1. To compare the accuracy of imaging markers measured by QyScore<sup>®</sup>, an FDA and CE marked medical device, and radiological assessment in distinguishing parkinson disease (PD) from progressive supranuclear palsy (PSP) patients
2. Any other indexes able to differentiate PD from PSP?



EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*

JULY 13-17, VIENNA & ONLINE

MYESR.ORG

# Methods

- 9 PSP patients
- 18 PD patients

- 25 healthy patients

52  
T1-3D brain MRIs



EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*  
JULY 13-17, VIENNA & ONLINE

MYESR.ORG

# Methods – radiological evaluation

M surface	3rdV	GcerbA
P surface	FH	BGA
M/P	MRPI	GCA
MCP	MRPI 2.0	MTA
SCP		



EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*  
JULY 13-17, VIENNA & ONLINE


MYESR.ORG





AR: 75,71 mm²  
MED: 739,01  
DS: 0,60  
Intervallo: 117,75 - 1062,79


AR: 511,52 mm²  
MED: 850,04  
DS: 62,91  
Intervallo: 33,21 - 1262,06



*BUILDING BRIDGES*  
JULY 13-17, VIENNA & ONLINE

MYESR.ORG






**NG BRIDGES**  
VIENNA & ONLINE


MYESR.ORG





**IDGES**  
& ONLINE

MYESR.ORG



# Methods – automated analysis

**1<sup>st</sup> step:** QyScore<sup>®</sup> software automatically segmented 17 brain structures, providing volumes and population-normed z-scores. The accuracy of the latter was compared with visual radiological assessment performed by an expert neuroradiologist;

**2<sup>nd</sup> step:** An automatic MRPI was developed by Qynapse and compared with the visual MRPI calculated by an expert neuroradiologist;

Metrics have been compared using Kruskal-Wallis test, Benjamini-Hochberg method. Overall diagnostic accuracy estimated as the area under the receiver operator curve with 95% CI




EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*



JULY 13-17, VIENNA & ONLINE

MYESR.ORG






*BUILDING BRIDGES*  
JULY 13-17, VIENNA & ONLINE

MYESR.ORG








# Results – 1<sup>o</sup> step

PD vs PSP	
Radiological evaluation	QyScore evaluation
Surface M (p=0.003)	Brainstem (p<0.001)
<b>M/P (p=0.005)</b>	Globus Pallidus (p<0.001)
MRPI (p=0.005)	Thalamus (p<0.001)
<b>MRPI2.0 (p=0.003)</b>	Amygdala (p<0.001)
PCM Diameter (p=0.046)	
PCS Diameter (p=0.033)	


	PD vs PSP	AUC
QyScore markers	Brainstem	0.9444
	Globus Pallidus	0.9603
	Thalamus	0.9524
	Composite (BS+GB+TH)	<b>0.9683</b>
Manual Indexes	Surface M	<b>0.8929</b>
	MRPI2.0	0.8889



p-value € **0.3085**



# Results – 1<sup>o</sup> step

	PD vs PSP	Sensitivity	Specificity
QyScore® markers	Brainstem z-score	0.857	0.944
	Globus Pallidus z-score	1.000	0.944
	Thalamus z-score	1.000	0.944
Radiological Assessment	MRPI2.0	1.000	0.722
	MRPI	0.714	0.888



- Globus Pallidus AUC = 0.96
- Brainstem AUC = 0.94
- Thalamus AUC = 0.95
- MRPI2.0 AUC = 0.88
- MRPI AUC = 0.77

# Results – 2<sup>o</sup> step




# Results – 2<sup>o</sup> step

## ACCURACY of QyScore markers and Visual Indexes

PD vs PSP			
	AUC	Sensitivity	Specificity
Manual MRPI	0.7698	0.7142857	0.8888888
Automatic MRPI	0.9286	0.8571429	0.9444444

— Manual MRPI  
— Automatic MRPI



# Conclusions

- Automated markers quantified using QyScore® as well as the automatic MRPI equally performed as an expert neuroradiologist in distinguishing PD and PSP patients.
- The radiological evaluation is a time consuming process, prone to the clinician's expertise and to inter-observer variability.
- AI and machine learning will allow to obtain precise and reproducible measures.



EUROPEAN CONGRESS OF RADIOLOGY

ECR 2022

*BUILDING BRIDGES*

JULY 13-17, VIENNA & ONLINE

MYESR.ORG

**Thank you for your attention**



**University of Verona**

Mambrin F.

Mansueto G.

Tinazzi M.

Di Vico I.

Pizzini F.B.



**QYNAPSE**

Cavedo E.

Guizard N.

Ayoub G.

EUROPEAN CONGRESS OF RADIOLOGY

**ECR 2022**

**BUILDING BRIDGES**

JULY 13-17, VIENNA & ONLINE

MYESR.ORG